



## Kit RAL 555

The kit RAL 555 is a quick variant of the Romanowsky stain. It is a multi-purpose stains that can be used in hematology applications (blood and bone marrow smears), parasitology, pathology, spermiology and cyto-bacteriology.

Most often used in hematology, this kit allows a reliable and reproducible staining of blood smears in 15 seconds.

It is also commonly used for the detection of malaria or the morphological study of spermatozoa (spermocytogram).


### KEY FEATURES

**Easy to use**

**Fast**

**Reproducible**


**Versatile**




555 staining of peripheral blood smear from patient with anormal CBC



555 staining - Spermogram



555 staining - Pleural fluid



555 staining - Malaria

## Kit RAL 555

### AVAILABLE KIT



Kit RAL 555

REF.: 361550-0000

### TECHNICAL SPECIFICATIONS

#### APPLICATION AREAS

- Hematology
- Parasitology
- Pathology
- Cyto-bacteriology
- Spermiology

#### NUMBER OF TESTS

- 100 to 200 slides

#### NUMBER OF BOTTLES

- 3 bottles

#### PRODUCT DESCRIPTION

- FIX RAL: Clear purple solution  
REF.: 362870 - 1 X 100 mL
- EOSIN RAL: Clear orange red solution  
REF.: 361640 - 1 X 100 mL
- BLUE RAL: Clear dark blue solution  
REF.: 361650 - 1 X 100 mL

#### REFILLS

- FIX RAL  
REF.: 362870-1000 - 1 X 1L  
REF.: 362870-2500 - 1 X 2.5L
- EOSIN RAL  
REF.: 361640-1000 - 1 X 1L  
REF.: 361640-2500 - 1 X 2.5 L
- BLUE RAL  
REF.: 361650-1000 - 1 X 1L  
REF.: 361650-2500 - 1 X 2.5 L