


## CellaVision® Remote Review Software

CellaVision Remote Review Software gives remote users access to processed slides and makes it possible to review and verify cell differentials from anywhere, anytime. By enabling remote database access, the software improves communication and collaboration between colleagues and makes it possible for morphology experts to access differentials remotely and provide consultation immediately.

When implemented together with CellaVision hardware and supporting applications, CellaVision Remote Review Software removes restrictive geographical constraints from the smear review process, enabling laboratories to work smarter and perform better.



### Features

Access cell images, both current and historical, remotely

Collaborate within / between sites and lab units in real time

View cell counter results and flags


## How does it work?

- Available as a single user license as well as concurrent user licenses
  - **Professional Edition:** for 1 single user
  - **Team Edition:** for up to 5 concurrent users
  - **Group Edition:** for up to 25 concurrent users
  - **Enterprise Edition:** for up to 75 concurrent users
- The single user license is installed directly on the users' PC
- The concurrent user licenses are installed on the hospital server or on multiple users' PC:s
- After installation, users can review, adjust and verify cell differentials remotely


## Why CellaVision Remote Review Software?

The Software creates workflow efficiency within and between labs by enabling remote access to processed slides.

- Removes geographical constraints from the smear review process
- Makes it possible to review and verify cell differentials from anywhere, anytime
- Facilitates collaboration between colleagues
- Facilitates review and consultation by morphology experts
- Creates an opportunity for flexible and cost-effective utilization of staff and resources



One WBC gallery view



Two WBC galleries view

## CellaVision® Remote Review Software

## Technical Specifications

### EDITIONS

- Professional Edition: for 1 single user
- Team Edition: for up to 5 concurrent users
- Group Edition: for up to 25 concurrent users
- Enterprise Edition: for up to 75 concurrent users

### SUPPORTED PRODUCTS

- CellaVision® DM1200
- CellaVision® DM9600
- CellaVision® DC-1

### SUPPORTED APPLICATIONS

- CellaVision® Peripheral Blood Application
- CellaVision® Advanced RBC Application
- CellaVision® Body Fluid Application

### SPECIFICATIONS

- Minimum monitor resolution: 1280 x 1024
- CellaVision® Remote Review is supported by the latest Windows® Operating systems.
- CellaVision® Remote Review is supported by the latest Citrix Virtual Apps and Desktops.
- Please refer to our IT Configuration Guidelines for further specifications.

Specifications are subject to change without notice. This product may not be available in all markets.